

user guide


LIPOVACUUM™



If you have any questions, or if you experience any difficulties please call our
Rio Careline on: **01242 702345** or email us at: **online@riobeauty.com**

LIPOVACUUM™ - SIMPLE YET EFFECTIVE

Congratulations on your purchase of Lipovacuum. Lipovacuum allows you to enjoy the benefits of Lipo-Dermology™ without the inconvenience of regular and expensive visits to a salon. Lipo-Dermology™ is a relatively simple treatment using vacuum created pulses which pass through the skin in waves, to repeatedly compress and stimulate the underlying tissue structures. This stimulating action helps improve skin elasticity, skin tone and reduce sagging skin, fine lines and wrinkles.



Used on body areas with localised subcutaneous fatty deposits (the skin often has an orange peel appearance) the pulsing compression waves work to break down the hardened fat cells resulting in a smoother and more toned appearance to the skin. The stimulation of blood supply and lymphatic drainage help bring an immediate supply of valuable nutrients and oxygen deeper into the affected tissues, while also carrying away released toxins.

The various treatment heads enable use on most areas of the body. Lipovacuum can be used on the face with the smallest head to target wrinkles and sagging skin around the jawline and neck.

Regular treatments are required for all areas to improve and maintain skin condition and tone.

If you have any questions please call our customer care representatives on **01242 702345**. They will be delighted to help you get the best results from your Lipovacuum.


Read all instructions thoroughly and retain for future reference.

GETTING TO KNOW YOUR UNIT

- 1 Lid
- 2 Hose Connector
- 3 Small Lift Treatment Head
- 4 Large Roller Treatment Head
- 5 Medium Roller Treatment Head
- 6 Filter
- 7 Hose
- 8 Hose Connection Point
- 9 Lid Catch
- 10 Mains Cable
- 11 On/Off Button
- 12 Power Increase Button
- 13 Power Indicator Lights
- 14 Power Decrease Button
- 15 Mode Button
- 16 Start, Stop and Pause Button
- 17 Filter Warning Light




TREATMENT HEADS




Large Roller Head

This head rolls over the skin during treatment. Use for large areas of the body, for instance, the legs, buttocks and back.



Medium Roller Head

This head also rolls over the skin during treatment and can be used on smaller areas of the body such as the arms.




Small Lift Head

This head uses a direct lifting action for precise work on small areas such as the face, neck and hands.

PRECAUTIONS


- Seek medical advice from your doctor before use if you are pregnant or you have any of the following: pacemaker, heart disorder, epilepsy, diabetes, multiple sclerosis, or if you are recovering from an operation.
- Do NOT use on scars, broken skin, rosacea, pustular or inflamed acne, abraded skin eczema, damaged, inflamed or diseased skin, or any other skin ailment. It is recommended that you consult a doctor if you are unsure.
- Do NOT use where capillary vessels can be seen or excessively on one area per treatment.
- Do NOT use within 12 hours of the following treatments: sauna, steam bath, sunbed, electrolysis, facial waxing.
- Do NOT use over piercings or jewellery.
- Do NOT use in the bath or shower. Do NOT immerse the unit or heads in water.
- This appliance is not intended for use by persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- NOT for the use of children under the age of 14. Keep out of reach of children and pets.
- Single person use only. Do NOT share the treatment heads. You can purchase spare treatment heads separately.
- Discontinue use in the event of irritation.
- For external use only.

INSTRUCTIONS FOR USE




STEP 1

Connect the hose to the hose connection point and the selected treatment head by simply pushing until it clicks into place.



STEP 2

Connect the mains cable to the back of the unit and into your power supply. Position your unit on a flat, firm surface.



STEP 3


Press the on/off button to turn the unit on.

The blue lights around the start button will flash.




STEP 4

When you are ready to begin press the start button, the light will stay on, and the vacuum pump will start. The pump will run for 60 minutes before automatically stopping.



STEP 5


The unit will default to the Fatty Tissue Breakdown mode. Press the mode button to cycle through the 4 modes and select the one you require.



STEP 6


Use the power + or – buttons to select the power level.

INSTRUCTIONS FOR USE




STEP 7

Move the treatment head along the skin in small steps moving the head between pulses. Keep a light pressure to achieve a vacuum seal. Do NOT press too hard.



DISCONNECTING THE HEADS

Press the release tab on the end of the hose and pull the connected treatment head off. Ensure that your fingers do not obstruct the tab movement of the connector on the underside of the hose.



STEP 8

Follow the advice in this user guide for the area of the body you are treating.



PAUSE

Press the start button to pause the unit at any time, the blue lights will flash until you press the button again to resume your treatment.



STEP 9

When you have finished turn off by pressing the on/off button. Store all the treatment heads and the hose back on the unit, then close the lid. Unplug the mains cable and store away from children.




POWER ON/OFF


The unit can be turned off at any time by pressing the on/off button.

NOTE: If no buttons are pressed after 20 minutes, the unit will enter energy saving mode. To restart, return to step 3.

CHANGING THE FILTER



The life expectancy of the filter is approximately 50 hours. If the red filter light is flashing the filter needs replacing. You can still use the Lipovacuum, but we recommend you change the filter as soon as possible. Contact Rio Customer Care on 01242 702345 or email online@riobeauty.com for replacement filters.



STEP 1

Plug the mains cable into the unit and into your power supply. Press the on/off button. The filter light will flash red.

Important: Keep the unit switched on during this procedure.




STEP 2

Twist the filter anti-clockwise until the outline circle lines up with the line on the unit. Pull the filter out of the unit.



STEP 3

Insert the new filter and twist clockwise until the solid circle lines up with the line to lock into place.




STEP 4

The filter warning light will now turn off.


Note: If the filter warning light is on, but not flashing, the filter is not fitted correctly. The unit will not work without a filter in place or if the filter is incorrectly fitted. In this case, remove and reinsert the filter correctly.

TREATMENT MODES




Fatty Tissue Breakdown

The pulsing action targets the subcutaneous layer, promoting the breakdown of the fatty tissue, commonly called cellulite. With regular sessions this can dramatically reduce the appearance of 'orange peel' skin.




Skin Toning and Firming

The extended pulsing action gently pulls the skin taut while stimulating blood circulation and providing a toning action.



Muscle Relief Massage

The massaging action of the vacuum gently works the muscles, leaving them feeling looser and relaxed.



Vascular Stimulation

The continuous suction stimulates the skin and increases blood flow leaving the skin invigorated and revitalised.

During Use


When using the Lipovacuum it is important to relax and take your time to work methodically over the area you are treating. For maximum benefit, try to maintain a vacuum seal with your skin while moving the treatment head over the skin in small steps inbetween pulses.

Follow the instructions in this user guide for the various modes and areas of the body. Whilst the Lipovacuum is a very safe and effective treatment it is important to remember the rules which apply to all massage products;

- Do NOT use over the same area repeatedly as this may result in skin soreness. It is important to also 'listen' to your own body feedback and if the treatment becomes uncomfortable reduce the power or move to another area of the body.
- Do NOT use over varicose veins, raised moles or growths.
- NEVER use directly on the spine, or back of knee and elbow joints. You should never apply pressure in these areas.
- When treating the eye area keep the treatment head at least 2 centimetres away from the eyes at all times. Do NOT place the treatment head immediately above or below the eye.

Finally, to aid lymphatic drainage always work the last movements in a sweeping direction towards the heart.

HIPS AND BUTTOCKS



The hips and buttocks tend to be one of the major areas of the body where fat is stored. Over time the build up of the fatty tissue deposits can develop to such a degree that the orange peel effect is visible on the surface of the skin. The breakdown of these hardened fat deposits are difficult to target

with diet and exercise alone. Lipovacuum can be used to directly target this layer of subcutaneous fat.

The lifting and rolling actions of Lipovacuum also stimulates the lymph system and local circulation. Encouraging lymphatic drainage helps to clear toxins and prevent the build-up of fatty deposits.



Only use on the areas shown, do not use directly on the fold under the buttock.



Fatty Tissue Breakdown - Large Head

Maximum 8 minutes on each side

Power level - up to maximum

Move in one direction keep the pressure firm but not painful, repeat over the area in the opposite direction. Use the highest power level which is comfortable to you.



Skin Toning and Firming - Large Head

8 minutes on each side

Power level - medium

Slowly and methodically work over the whole area in one direction only, strokes should flow towards the heart.



Muscle Relief Massage - Large Head

Maximum 8 minutes on each side

Power level - comfortable to maximum

Work the side and large area of the buttock in a left to right motion, followed by an up and down motion.





Vascular Stimulation - Large Head

8 minutes on each side

Power level - comfortable

Work the side and large area of the buttock in a left to right motion, followed by an up and down motion. Use a power level which is comfortable to you.

THIGHS



The thighs, along with the buttocks, are the one of the most common areas where the body naturally stores fat.

As with the hips and buttocks, the lifting and rolling actions of the Lipovacuum promotes lymphatic drainage combined with the massaging action to help prevent the build-up of toxins and fatty deposits whilst encouraging the breakdown of existing subcutaneous fat.

With the thighs containing some of the largest muscle groups in the body, muscle aches and pains are common. The Muscle Relief Massage programme is ideal for treating these robust muscles and providing rapid relief from tired, aching muscles.



Never use the Lipovacuum over the back of the knee. Do NOT use over areas of varicose veins.



Fatty Tissue Breakdown - Large/Medium Head
Maximum 8 minutes on each side
Power level - up to maximum

Move in vertical strokes, keeping the pressure firm, repeat over the area in lateral strokes. Finish with vertical strokes towards the heart. Use the highest power level which is comfortable to you.



Skin Toning and Firming - Large/Medium Head
6 minutes on each side
Power level - medium

Slowly and methodically work over the whole area with upward vertical strokes towards the heart.



Muscle Relief Massage - Large/Medium Head
Maximum 8 minutes on each side
Power level - comfortable to maximum


Move in vertical strokes, keeping the pressure firm, repeat over the area in lateral strokes. Finish with vertical strokes towards the heart.



Vascular Stimulation - Large/Medium Head
6 minutes on each side
Power level - comfortable

Move in vertical strokes, keeping the pressure firm, repeat over the area in lateral strokes. Finish with vertical strokes towards the heart. The power level should allow easy movements.

ABDOMEN



The abdomen has no bone structure protecting the major organs. It is, therefore, very important to only use gentle pressure and a low power on this area.

There are three layers of abdominal muscles: external oblique, internal oblique and the transversus abdominis, supported on the sides by rectus abdominis. These are regarded as the muscles responsible for core stability - vital for good posture and defining our shape.

Fat deposits around the abdominal area often build up slowly over time and, consequently, respond slowly to diet and exercise alone. Using the Lipovacuum helps to speed up the breakdown of the fatty tissue enhancing the effects of diet and exercise.

If you suffer from digestive discomfort you can use the Vascular Stimulation programme and work in a clockwise direction working from the outer area inwards. This follows the passage of the intestines and can help relieve some discomfort.



Do not use over the belly button area of the stomach. If in doubt consult your doctor for advice.



Fatty Tissue Breakdown - Large Head

Maximum 5 minutes

Power level - low

Move horizontally across the area with very light pressure. Repeat with vertical strokes. Keep the power level low and do not exceed the 5 minute time limit.



Skin Toning and Firming - Large/Medium Head

Maximum 5 minutes

Power level - low

Move horizontally across the area with very light pressure. Repeat with vertical strokes. Keep the power level low to allow easy movements.



Muscle Relief Massage

We do not recommend using this programme on the abdomen.




Vascular Stimulation - Medium Head

Maximum 5 minutes

Power level - low to comfortable

Ideally use a clockwise circular movement over the area. Use a low power level which is comfortable to you.

BACK



The back contains three major muscle groups: the latissimus dorsi, located on each side of your back, the erector spinae (lower back) which is made up of three muscles that run the full length of your back, and the rhomboids between the shoulder blades. These muscles are vital to our posture, mobility and our very quality of life. The back can easily become tense and knotted contributing to, and even causing, headaches, neck aches, bad posture, back pain and shallow breathing.

You will need a friend to use the Lipovacuum on your back. Lie face down on a comfortable surface and concentrate on relaxing and breathing in a rhythmic way.

We recommend using the large roller head, although some people may find the medium head more comfortable. Start at the shoulders, working on one then the other. Move to the back and use sweeping strokes from the centre outwards. Finish with long strokes up the back towards the head.



Do NOT use directly on the spine. Do NOT use over raised moles.



Fatty Tissue Breakdown - Large Head

8 minutes on each side

Power level - up to maximum

Move slowly keeping the pressure firm but not painful. Use the highest power level which is comfortable to you.



Skin Toning and Firming - Large Head

8 minutes on each side

Power level - medium

Slowly and methodically work over the whole area with a medium power level.



Muscle Relief Massage - Large Head

Maximum 8 minutes on each side

Power level - comfortable to maximum

Move slowly keeping the pressure firm but not painful. Use a power level which is comfortable to you.




Vascular Stimulation - Large Head

Maximum 10 minutes on each side

Power level - comfortable

Move slowly keeping the pressure firm but not painful. Use a power level which is comfortable to you.

SHOULDERS



The shoulder muscles can easily become tense and knotted, this tension can contribute to headaches, neck aches, bad posture and inhibit relaxing breathing which in turn can increase stress. Using the Lipovacuum over the shoulder area relaxes the muscles.

During use try to concentrate on dropping the shoulders and breathing in a rhythmic way to fully benefit from the action. Move the treatment head back and forth over the shoulder trying to maintain contact with the skin at all times.

Although you can use the Lipovacuum on your own shoulders, one at a time, having a friend to help will allow you to relax more and feel the full benefits.



Do NOT use directly on the spine or back of the neck.
Do NOT use over raised moles.



Fatty Tissue Breakdown - Medium Head
4 minutes on each side
Power level - up to maximum

Move across the shoulder in one direction keeping the pressure firm but not painful, repeat over the area in the opposite direction. Use the highest comfortable power level.



Skin Toning and Firming - Medium Head
Approximately 4 minutes on each side
Power level - medium

Slowly and methodically work over the whole area in one direction covering all the skin.



Muscle Relief Massage - Medium Head
Maximum 6 minutes on each side
Power level - comfortable to maximum


To relieve the tension in shoulder muscles move the medium head across the shoulder in both directions. Use a power level which is comfortable to you.



Vascular Stimulation - Medium Head
4 minutes on each side
Power level - comfortable

Move the medium head across each shoulder in both directions. Use a power level which is comfortable to you.

DECOLLAG



Just like our faces, the decolage (the area above the breasts to the neck) is an area of thin skin which shows the stress and strain of daily life very easily. The skin texture becomes less toned over time eventually giving a 'crepe' appearance. Using the

Skin Toning and Firming programme on a low power and with light pressure provides a gentle relaxing treatment stimulating blood circulation through the skin improving the complexion, toning, firming and reducing fine lines and wrinkles.

It is important not to 'over-work' this area as the skin is delicate and can easily bruise.

Do not finish with strokes towards the heart for this area of the body.



Do not use if you have a history of cardiac or chest problems. Consult your doctor for advice.



Fatty Tissue Breakdown

We do not recommend using this programme on the decolage.



Skin Toning and Firming - Small Head

12 minutes

Power level - comfortable

Slowly and methodically work over the whole area following the directions in the illustrations. The power level should allow easy movements.



Muscle Relief Massage

We do not recommend using this programme on the decolage.




Vascular Stimulation - Medium/Small Head

12 minutes

Power level - comfortable

Slowly and methodically work over the whole area following the directions in the illustrations. The power level should allow easy flowing movements.

CALF MUSCLES



The calf muscles are often a neglected area of the body, yet they take a huge amount of stress and strain everyday. This is especially true for people whose day involves a lot of standing.

Twinges and aches are often ignored until they develop into problems. Using the Lipovacuum regularly, muscle tensions can be treated at an early stage and serious problems can be avoided.

All the Lipovacuum programmes promote blood flow to the skin surface increasing circulation. Whilst this is beneficial all over the body, this is particularly good for the calves. (Do NOT use where varicose veins are already present.)

Some fitness programmes, for example, running or long distance cycling, can result in 'tender' calf muscles. Whilst the Vascular Stimulation programme can help alleviate the discomfort, it is important to use a low power level and treat the muscle very gently.



Never use the Lipovacuum over the back of the knee. Do NOT use over areas of varicose veins.



Fatty Tissue Breakdown - Medium Head

4 minutes on each calf

Power level - comfortable

Use vertical strokes over the muscle area keep the pressure firm but not painful, repeat with lateral strokes. Finish with more vertical strokes towards the heart.



Skin Toning and Firming - Medium Head

4 minutes on each calf

Power level - comfortable

Use vertical strokes towards the heart. The power level should allow easy movements.



Muscle Relief Massage - Medium Head

4 minutes on each calf

Power level - comfortable

Use vertical strokes over the muscle area, repeat with lateral strokes. Finish with more vertical strokes towards the heart.





Vascular Stimulation - Medium Head

4 minutes on each calf

Power level - comfortable

Use vertical strokes towards the heart. The power level should allow easy movements.

ARMS



Our arms and hands are used constantly throughout the day resulting in tension and sometimes pain. The Lipovacuum can massage this tension

away improving our sense of well-being.

The upper arm area is also a source of concern to many women who suffer from 'bingo-wings'. This is often a combination of stored fat and flaccid muscle which sags from the upper arm. Using the Lipovacuum regularly on this area will help to firm and tone up the skin on this area.



Remove all items of jewellery on the arms before using the Lipovacuum in this area.



Fatty Tissue Breakdown - Medium Head
4 minutes on each arm
Power level - up to maximum

Use vertical strokes over the lower then upper arm area, keeping the pressure firm but not painful, repeat with lateral strokes. Finish with long strokes towards the heart.



Skin Toning and Firming - Medium Head
4 minutes on each arm
Power level - medium

Slowly and methodically work over the whole area with vertical strokes towards the heart.



Muscle Relief Massage - Medium Head
4 minutes on each arm
Power level - comfortable to maximum


Use vertical strokes over the lower then upper arm area, keeping the pressure firm but not painful, repeat with lateral strokes. Finish with long strokes towards the heart.



Vascular Stimulation - Medium Head
4 minutes on each arm
Power level - comfortable

Use vertical strokes towards the heart. The power level should allow easy movements.

FACE AND NECK



Our faces are one of the areas which quickly shows the stress and strain of daily life. Lines and wrinkles become more pronounced as we adopt the same expressions day after day. Using the Lipovacuum on a low power and with very light pressure provides a gentle relaxing treatment which helps to release toxins, stress and tension while plumping up the skin to reduce the appearance of fine lines and wrinkles.

As we grow older the facial muscles lose tone and the skin loses its elasticity. As a result, wrinkles, laughter lines and crow's feet develop and the skin can look dull and lifeless. The lifting action of the Lipovacuum stimulates blood circulation to the skin bringing oxygen to the tissues improving the complexion, toning, firming and reducing fine lines and wrinkles, to rejuvenate your skin.

By using the small lifting treatment head on an individual wrinkle Lipovacuum can be used to target specific areas.



Always clean your skin before use. Keep the treatment head at least 2 centimetres away from the eyes at all times. Do NOT place the treatment head immediately above or below the eye.

LIPOVACUUM™



Fatty Tissue Breakdown

We do not recommend using this programme on the face and neck.



Skin Toning and Firming - Small Head

10 minutes

Power level - low to comfortable

Slowly and methodically work over the whole area following the directions in the illustrations. The power level should remain very low. Do not 'drag' the skin.



Muscle Relief Massage

We do not recommend using this programme on the face and neck.



Vascular Stimulation - Small Head

10 minutes

Power level - comfortable

Slowly and methodically work over the whole area following the directions in the illustrations. The power level should remain very low to allow easy movements.

HANDS



The hands and fingers are used constantly throughout the day. Massage in this area can help to release some of the built-up tension and help to relax the muscles.

The massaging and lifting actions of the Lipovacuum help to improve blood circulation to the hands and fingers.

The massaging actions give relief from aches and stiffness and can help to prevent conditions such as RSI (Repetitive Strain Injury).



Remove all items of jewellery on the arms before using the Lipovacuum in this area.



Fatty Tissue Breakdown

We do not recommend using this programme on the hands.



Skin Toning and Firming - Small Head

Up to 5 minutes on each hand

Power level - low to comfortable

Slowly and methodically work over the whole area following the directions in the illustrations. Finish with strokes from the finger tips towards the wrist.



Muscle Relief Massage - Small Head

Up to 5 minutes on each hand

Power level - comfortable

To relieve aching hands work over the palm area (the soft muscle area) with a light pressure in both directions. Use a power level which is comfortable to you.





Vascular Stimulation - Small Head

Up to 5 minutes on each hand

Power level - comfortable


Slowly and methodically work over the whole area following the directions in the illustrations. Finish with strokes from the finger tips towards the wrist.

STORING THE HOSE




STEP 1

Place one end of the hose in the hose recess at the top right corner of the unit.



STEP 2

Push the hose into the recess around the outer edge of the unit.



STEP 3

The other end of the hose will finish at the top left hand corner of the unit. As you close the lid, ensure the hose end is tucked inside.

WARNINGS AND MAINTENANCE

- Only use the product with the mains lead supplied.
- Do NOT leave the unit unattended when plugged in.
- If the mains cable is damaged do NOT use it. In the event of a fault a suitable replacement should be used.
- After use, unplug the unit and store out of reach of children.
- There are no user serviceable parts inside the unit.
- Misuse of the unit can cause harm.
- If you experience unexpected behaviour of this product then it is possible some external electromagnetic interference is responsible. You should remove this source and restart the product to restore full operation.
- Should you notice a brief white flash at the power level indicator lights, this is part of the internal self-check and can be ignored.
- Do NOT place the unit near any heat sources.
- Do NOT use outdoors. Do NOT use the unit where aerosols or spray products are being used. Do NOT use the unit where oxygen is being administered.
- Store the unit at room temperature. Do not store at temperatures below 5°C.
- Do NOT carry the unit by the lid or mains cable.
- Coil the hose back neatly after each use.
- To clean, unplug the unit and use a slightly damp sponge to wipe the outside and inside under the lid. Do NOT use solvents or detergents.
- The treatment heads may be cleaned by wiping with a damp cloth or antiseptic wipe. Do NOT use solvents or detergents.



THE UNIT CONFORMS WITH CURRENT REGULATORY STANDARDS FOR SAFETY AND ELECTROMAGNETIC COMPATIBILITY



INDOOR USE ONLY

CUSTOMER CARE

If you have any questions, or if you experience any difficulties please call our Rio Careline on: **01242 702345**
or email us at: **online@riobeauty.com**

GUARANTEE

To register your two year guarantee online, visit: www.riobeauty.com/warranty_online.htm

DISPOSAL



Electrical and electronic waste is the fastest growing waste stream in the UK. Around 1.8 million tonnes are generated every year. The Waste of Electrical and Electronic Equipment (WEEE) Directive (2002/96/EC) has been put in place to recycle products using best available recovery and recycling techniques to minimise the impact on the environment, treat any hazardous substances and avoid the increasing landfill.

The product is classed as Electrical or Electronic equipment so please ensure that at the end of the life of this product it is disposed of in the correct manner in accordance with local authority requirements, or, returned to The Dezac Group Ltd., Cheltenham, PO Box 54, GL14 2WW, UK. It must NOT be disposed of with household waste.

Rio health, beauty and fitness products incorporate the latest technology to help you achieve the best results.

To ensure complete satisfaction we constantly review the performance of our products working in partnership with customers and researchers.

As part of this programme we welcome your comments and suggestions which can be sent to us at the address below.

www.riobeauty.com

UK



Rio, Dezac, Lipovacuum and Lipo-Dermology are trademarks of The Dezac Group Ltd
© The Dezac Group Ltd 2010. Manufactured by The Dezac Group Ltd, PO Box 17, Cheltenham, England GL53 7ET

1-IN-CEVC-UK/2.1